INDICATIONS FOR USE:


The T-Piece Resuscitator Circuit is a component used in conjunction with T-Piece resuscitator devices.

ACCESSORIES:

Silicone Contour Masks Sizes: #0 and #1
Cushion Mask Sizes: Neonatal and Infant

WARNINGS:

1. Incorrect operation of this circuit without a manometer and pressure relief can be hazardous.
2. Do not use in the presence of flammable materials and ensure that no sources of ignition are present while in use. Fire hazards are possible in oxygen enriched environments.
3. Users should release the T-piece (manual control) as soon as the breath is delivered. Failure to do so will result in an extended breath, which may prevent exhalation.
4. Do not leave patients unattended when in use.
5. The recommended gas flow range is from 5 – 15 LPM. DO NOT EXCEED 15 LPM.
6. Remove protective cap from the T-Piece circuit prior to connecting the T-Piece to the mask and before patient use (see **Figure 1 & Figure 2**).

**CAUTIONS:**

1. Never wait to begin mouth-to-mask resuscitation if a resuscitator is not immediately available or cannot be used effectively (see your department manual or AHA guidelines for accepted patient resuscitation procedure).
2. US Federal Law restricts this device to sale or on the order of a physician.
3. To be used only by persons who have been trained in infant/neonatal resuscitation.
4. Any T-Piece resuscitator system should be used only after completion of a pre-use verification check to ensure that correct ventilation pressures will be delivered to the patient. Discard any circuit that fails a pre-use check.
5. Oxygen concentration should be monitored at all times using an oxygen analyzer.
6. Adjusting input flow rate will affect PIP and PEEP.
7. Do not attempt to disinfect any part of this circuit.
8. If patient valve becomes occluded shake or rinse with water.

PRE-USE CHECK:

1. Refer to the manufacturer of the resuscitator control units operating instruction for pre use check.

DIRECTIONS FOR USE:

1. Refer to the manufacturer of the resuscitator control units operating instruction for directions for use.

SPECIFICATIONS:

1. Dead space: 4 ml
2. Expiratory Resistance: 0.2 cm H₂O at minimum PEEP setting @ 6 LPM
3. Inspiratory Resistance: 1.6 cm H₂O at minimum PEEP setting @ 6 LPM
4. Storage Temperature Range: - 40°C to + 60°C, up to 95% humidity
5. Operating Temperature Range: - 18°C to + 50°C, up to 95% humidity
6. Patient connection ISO 5356-1 – 15 mm taper female
7. Gas Delivery Pressures:
Positive End Expiratory Pressure (PEEP) range at:
@ 5 LPM up to approximately 2 cm H₂O
@ 8 LPM up to approximately 6 cm H₂O
@ 10 LPM up to approximately 9 cm H₂O
@ 15 LPM up to approximately 15 cm H₂O

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*Registered trademark in USA for Mercury Enterprises, Inc.