


Available Now . . .

New CPR-2 Resuscitation System



Finally,

a more affordable resuscitation system with the Mercury Medical quality you can count on.

MercuryMedical

Mercury Medical's CPR-2 Bags with Optional Color-Coded Manometer

The CPR-2 is a full line of resuscitators manufactured with new thermoplastic bag material providing superior lung compliance “feel” with exceptional responsiveness. The CPR-2 now offers clinicians a choice of the same unparalleled quality performance and accessories as the regular synthetic rubber CPR bag product line.



Squeeze the resuscitator bag to ventilate and the manometer verifies airway pressure.



Release squeeze on exhalation and the manometer verifies PEEP pressure (when PEEP valve is attached to the resuscitator).

BENEFITS	MERCURY MEDICAL	OTHER SYSTEMS
• Optional Color-Coded Manometer, accuracy is ± 3 cm H ₂ O between 0 - 15 cm H ₂ O	X	-
• Monitors both Airway and PEEP Pressure	X	-
• Fast re-expansion/high breath-per-minute capability	X	-
• Latex-Free medical grade thermoplastic bag material	X	X
• Optional CO ₂ detector accessory	X	-
• Patient Valve detaches to add Omni-Link tubing that provides greater reach capability, ideal for neonatal and MRI needs	X	-
• Optional Exhalation Filter and MR Conditional Peep Valve accessory	X	-

The Resuscitation Bag with Color-Coded Manometer allows for in-line monitoring of both airway and PEEP pressures.

The Mercury Medical Disposable Color-Coded Manometer is ready when you need it.

Manufactured by:

Mercury Medical RESPIRATORY

Patent # US 5,557,049 B1

ISO 13485

Clearwater, FL 33762-4807 · 800.237.6418 · www.mercurymed.com

Copyright 02/09

ML2009-122

Rev. A

Mercury Medical's CPR-2 Bags with Optional Color-Coded Manometer



Ordering Information

Adult		
#10-56301	Adult CPR, No Mask, w/Manometer, Expandable Reservoir	12/Box
#10-56302	Adult CPR, w/Mask, w/Manometer, Expandable Reservoir	12/Box
#10-56303	Adult CPR, no Mask, Expandable Reservoir	12/Box
#10-56304	Adult CPR, w/Mask, Expandable Reservoir	12/Box
#10-56001	Adult CPR, No Mask, w/Manometer, O ₂ Reservoir Bag	6/Box
#10-56002	Adult CPR, w/Mask, w/Manometer, O ₂ Reservoir Bag	6/Box
#10-56003	Adult CPR, No Mask, O ₂ Reservoir Bag	6/Box
#10-56004	Adult CPR, w/Mask, O ₂ Reservoir Bag	6/Box
Child		
#10-56101	Child CPR, No Mask, w/Manometer, O ₂ Reservoir Bag	6/Box
#10-56102	Child CPR, w/Mask, w/Manometer, O ₂ Reservoir Bag	6/Box
#10-56103	Child CPR, No Mask, O ₂ Reservoir Bag	6/Box
#10-56104	Child CPR, w/Mask, O ₂ Reservoir Bag	6/Box
Infant		
#10-56201	Infant CPR, No Mask, w/Manometer, O ₂ Reservoir Bag	6/Box
#10-56202	Infant CPR, w/Mask, w/Manometer, O ₂ Reservoir Bag	6/Box
#10-56203	Infant CPR, No Mask, O ₂ Reservoir Bag	6/Box
#10-56204	Infant CPR, w/Mask, O ₂ Reservoir Bag	6/Box

CPR-2 Accessories

#10-55350	Single-Patient-Use Pressure Manometer 0 - 60cm H ₂ O	20/Box
#10-55356	Manometer with Tee, Left Orientation 0 - 60cm H ₂ O	20/Box
#10-55357	Manometer with Tee, Right Orientation, 0 - 60cm H ₂ O	20/Box
#10-55375	Exhalation Filter, 30mm O.D. x 30mm I.D.	20/Box
#10-55360	OmniLink [®] 12" Length	50/Box
#10-55361	OmniLink II, 6" Length	25/Box
#10-55330	PEEP Valve with 22mm Adapter (MR Conditional)	12/Box
#10-55370	StatCO ₂ [®] End Tidal CO ₂ Detector (for patients with body weight greater than 15 kg.)	20/Box
#10-55371	Mini StatCO ₂ [®] End Tidal CO ₂ Detector (for patients with body weight between 1 and 15 kg)	20/Box
#10-55373	Neo-Stat CO ₂ <Kg [™] For patients below 2 kg (0.25 kg - 2.0 kg)	20/Box



Color-Coded
Manometer



Exhalation
Filter



Omni-
Links



PEEP
Valve
MR Conditional



StatCO₂



Mini
StatCO₂



Neo-StatCO₂ <Kg