

ONLY ONE

Introducing a First in **PREEMIE** Face Mask Technology...



m a s k s

Mercury's industry's first mask design was developed with your needs in mind and offer the following advantages:

- Soft & flexible
- Ergonomic design
- Tight seal & reduced mask leak vs. cushion style masks
- Two styles to meet your preemie patient population

Larger mask sizes are also available, see inside for more details.



www.mercurymed.com

Mercury Medical®

Your Need . . . Our Innovation®

Introducing an **ONLY ONE** industry first in face mask technology for your preemie patients! Mercury's new anatomical silicone face masks with flange are available in sizes 43mm (Preemie) and 37mm (Ultra Preemie). So, instead of using the smallest masks that still cover the baby's eyes, or having to intubate when the smallest sizes are not available... Finally, the masks you've been waiting for! Mercury Medical's...

PREEMIE MASKS



New! Preemie Anatomical Silicone Mask with Flange (43mm)



Neonatal Cushion Mask Covering Eyes



Two equivalent tiny round silicone masks are also available so that you can have the best combo to meet the needs of any facial structure, whether the neonates have a flatter face or those that have a pronounced nasal bridge.



ORDERING INFORMATION		10/Box
Part #	Mask Description	Suggested Weight Range
10-50610	Anatomical Silicone 43mm Preemie (Very Low Birth Weight)	<1.5kg (3.3 lb)
10-50611	Anatomical Silicone 37mm Ultra Preemie (Extremely Low Birth Weight)	400g - 1kg (.88 - 2.2 lb)
10-50609	Round Silicone 42mm ('00') Preemie (Very Low Birth Weight)	<1.5kg (3.3 lb)
10-50612	Round Silicone 35mm ('000') Ultra Preemie (Extremely Low Birth Weight)	400g - 1kg (.88 - 2.2 lb)

In addition to our new **PREEMIE** masks, check out our other masks of various sizes, also featured in this brochure: the Disposable Anatomical Face Masks, the Reusable Anatomical Silicone Face Masks and the Reusable Round Silicone Face Masks. ➤

Disposable Anatomical Face Masks

The Disposable Anatomical Face Masks feature a cushion that's anatomically correct and contours for all patient sizes. They provide an overall tight seal and are designed specifically for all patient sizes... neonatal through adult. The masks are DEHP-FREE.

FEATURES	BENEFITS
Completely Disposable	Eliminates Cross-Contamination Risks and Associated Reprocessing Costs
Clear Material Construction	Better Patient Visualization
Anatomically Designed	Accommodates a Wide Range of Facial Features and Reduces Patient Discomfort
Soft, Pliable, Air-Filled Cushion	Reduces Wrinkles or Creases and Allows for a Secure, Low Pressure Seal
Removable, Color-Coded Hook Ring Fastener	Allows for Use of Traditional Headstraps



ORDERING INFORMATION

Part #	Description	Packaged
#10-50500	Adult Large Clear Anatomical Face Mask w/ 22mm Female Connector	30/Box
#10-50501	Adult Small Clear Anatomical Face Mask w/ 22mm Female Connector	30/Box
#10-50502	Child Clear Anatomical Face Mask w/ 22mm Female Connector	20/Box
#10-50503	Pediatric Clear Anatomical Face Mask w/ 22mm Female Connector	20/Box
#10-50504*	Infant Clear Anatomical Face Mask w/ 15mm Male Connector	20/Box
#10-50505*	Neonatal Clear Anatomical Face Mask w/ 15mm Male Connector	20/Box

*NOTE: Hook Ring Fastener not available on these sizes.

www.mercurymed.com

Reusable Anatomical Silicone Face Masks

Mercury Medical's Reusable Anatomical Silicone Face Masks are anatomically contoured to conform to natural facial profiles allowing for a secure, low pressure seal.



**Child
Size 3**



**Small Adult
Size 4**



**Adult
Size 5**

ORDERING INFORMATION		
Part #	Description	Packaged
#10-50605	Child Mask Size #3	10/Box
#10-50607	Small Adult Mask Size #4	10/Box
#10-50604	Adult Mask Size #5	10/Box

Reusable Round Silicone Face Masks



Size 0



Size 1



Size 2

ORDERING INFORMATION		
Part #	Description	Packaged
#10-50601	Round Type Silicone Face Mask, Size 0 (48mm) 15mm Male Connector	10/Box
#10-50602	Round Type Silicone Face Mask, Size 1 (60mm) 15mm Male Connector	10/Box
#10-50603	Round Type Silicone Face Mask, Size 2 (76mm) 22mm Female Connector	10/Box
#10-50608	Round Type Silicone Face Mask, Size 2 (76mm) 15mm Male Connector	10/Box

MercuryMedical[®]

11300 - 49th Street North • Clearwater, Florida 33762-4807 USA

ORDERING INFORMATION: 800.835.6633 or 727.573.0088

Fax: 727.571.3922

www.mercurymed.com