

ONLY ONE



NIF (to -40 cm H₂O)
10-340-00
20/Box

Single-Patient-Use

NIFometer®

NEGATIVE INSPIRATORY FORCE METER

The quicker, easier, safer, more reliable and convenient way to measure NIF.

INTRODUCING

NIFometer (to -60 cm H₂O)



10-340-02
20/Box

QUICK AND EASY

- Small, compact, self-contained. Connects directly to ET tube eliminating need for separate tubing and connectors.
- Includes built-in occlusion button that you simply depress during test.

SAFE

- Single-patient-use eliminates risk of cross-contamination.
- Provides unobstructed view of patient during test.

RELIABLE

- Accurate to ± 2 cm H₂O -20 to -30 cm H₂O (NIF 0-40) & ± 3 up to -15 cm H₂O & ± 5 from -15 to -60 cm H₂O (NIF 0-60).
- No calibration required.
- Includes memory indicator pointer that records highest force reached and can then be reset.

CONVENIENT

- Lightweight weighing only 30 Grams.
- Small and compact, fits in pocket – readily available for added efficiency.
- Also available in kit with elbow, mouthpiece and nose clip.

www.mercurymed.com

MercuryMedical®

Single-Patient-Use

NIFometer®

NEGATIVE INSPIRATORY FORCE METER

The quicker, easier, safer, more reliable and convenient way to measure NIF.



NIFometer
(to -40 cm H₂O)

NIFometer
(to -60 cm H₂O)



NIFometer Kit
Available with NIFometer
to -40 cm H₂O & to -60 cm H₂O

SPECIFICATIONS

Range	To -40 cm H ₂ O
Variance	± 2 cm H ₂ O -20 to -30 cm H ₂ O
Fitting	15 mm I.D./22 mm O.D.
Weight	30 grams nominal

Range	To -60 cm H ₂ O
Variance	± 3 up to -15 cm H ₂ O & ± 5 from -15 to -60 cm H ₂ O
Fitting	15 mm I.D./22 mm O.D.
Weight	30 grams nominal

ORDERING INFORMATION

Item	Description	Packaged
10-340-00	NIFometer (to -40 cm H ₂ O) Single-Patient-Use	20/Box
10-340-01	NIFometer Kit, (to -40 cm H ₂ O) Includes: elbow, mouthpiece & nose clip	10/Box
10-340-02	NIFometer (to -60 cm H ₂ O) Single-Patient-Use	20/Box
10-340-03	NIFometer Kit, (to -60 cm H ₂ O) Includes: elbow, mouthpiece & nose clip	10/Box

For more information about the NIFometer please call
toll-free **800.237.6418** or visit our website.