

Distributed by:  
**Mercury Medical**<sup>®</sup>  
*Your Need ... Our Innovation*<sup>®</sup>

# IGuard Eye Protector

## *Features:*

- Non-irritating
- Crystal clear lenses
- Easy tear-off peel backing
- Soft self-adhering foam



ARE YOUR PATIENTS' EYES  
**PROTECTED?**

**MercuryMedical**<sup>®</sup>

*Your Need ... Our Innovation*<sup>®</sup>

11300 - 49th Street North  
Clearwater, FL USA 33762-4807

PRSRT STD  
U.S. POSTAGE  
**PAID**  
PERMIT NO. 579  
ST. PETERSBURG, FL



Item #	Packaged
1091100	25/Box

**CALL TODAY**  
**1-800-835-6633**  
or visit us at  
[www.mercurymed.com](http://www.mercurymed.com)

# Face Masks

## Anatomical Face Masks



Mercury Medical's Face Masks are anatomically contoured to conform to natural facial profiles allowing a secure, low pressure seal without sharp edges that could damage the eyes.

### FEATURES:

- Completely disposable
- Clear material construction
- Anatomically designed
- Soft, pliable, air-filled cushion
- Removable color-coded headstrap fastener

### BENEFITS:

- Eliminates cross-contamination risks and associated reprocessing costs.
- Better patient visualization.
- Accommodates a wide range of facial features and reduces patient discomfort.
- Reduces wrinkles or creases and allows for a secure low pressure seal.
- Allows for use of traditional headstraps.

### ITEM NUMBER:

### DESCRIPTION:

10-50500	Adult Large Clear Anatomical Face Mask 30/case (Large)
10-50501	Adult Small Clear Anatomical Face Mask 30/case (Medium)
10-50502	Child Clear Anatomical Face Mask 20/case (Small)
10-50503	Pediatric Clear Anatomical Face Mask 20/case
10-50504*	Infant Clear Anatomical Face Mask 20/case
10-50505*	Neonatal Clear Anatomical Face Mask 20/case

\*Headstrap Fastener not available on these sizes.

## Round Silicone Masks

Features anatomically correct "average" contour for infant faces. Provides overall tight seal with internal cavity, minimizing dead space. Designed specifically for children under the age of ten (10).

10-50601	Size 0 10/pkg.
10-50602	Size 1 10/pkg.
10-50603	Size 2 10/pkg.



## PEEP VALVE

(MR Conditional)

The PEEP Valve works by applying an integral gas pressure to hold the alveolus open throughout inspiration and expiration. This internal pressure aids in the prevention of atelectasis and associated shunt conditions.<sup>1</sup>

### Features

- Built-in damping mechanism
- 30mm female fitting
- 0-20 cm H<sub>2</sub>O (mbar) adjustable range
- Single-patient-use
- Includes 22mm male adapter
- When using PEEP, the patient's condition should be monitored and PEEP pressure verified with a manometer.
- Use an MR system with a static magnetic field of 1.5 Tesla or less
- Use an MR system with a spatial gradient of 450 gauss/cm or less

### Benefits

- Reduces noise
- Conforms to ISO inlet standards

### ORDERING INFORMATION

Part #	Description	Packaged
#10-55330	PEEP Valve (disposable)	12/Box



1. Reference: Roger C. Bone, David H. Eubanks, Comprehensive Respiratory Care. The C.V. Mosby Co., St Louis. 1990;23:684.