

Mercury Medical's Latex-Free CPR Bags

Mercury Medical's CPR Bags are designed to accommodate the most urgent and unique situations of on-site EMS, Hospital I.C.U., E.R., O.R. and Recovery Room needs.



Adult CPR Bag

Patients Over 30 kgs

Maximum Rate	92 BPM		
Maximum Volume	900 - 1500 ml		
Pressure Relief Valve	(Optional)		
Oxygen Concentration			
LPM	3	5	10
Rate	12	12	12
Tidal Volume	500 ml	500 ml	500 ml
W/Reservoir	68 %	95 %	100 %
W/O Reservoir	34 %	47 %	70 %
Measured Volume			
Bag	1600 ml		
Reservoir	2600 ml		



Child CPR Bag

Patients Between 7 and 30 kgs

Maximum Rate	168 BPM
Maximum Volume	360 ml
Pressure Relief Valve	40 cm H ₂ O
Oxygen Concentration	
LPM	10
Rate	20
Tidal Volume	250 ml
W/Reservoir	100%
W/O Reservoir	85%
Measured Volume	
Bag	500 ml
Reservoir	2600 ml



Small Adult CPR Bag

Patients Over 30 kgs

Maximum Rate	110 BPM		
Maximum Volume	500 - 750 ml		
Pressure Relief Valve	(Optional)		
Oxygen Concentration			
LPM	3	5	10
Rate	12	12	12
Tidal Volume	500 ml	500 ml	500 ml
W/Reservoir	68 %	95 %	100 %
W/O Reservoir	34 %	47 %	70 %
Measured Volume			
Bag	1000 ml		
Reservoir	2600 ml		



Infant CPR Bag

Patients Under 7 kgs

Maximum Rate	180+ BPM
Maximum Volume	175 ml
Pressure Relief Valve	40 cm H ₂ O
Oxygen Concentration	
LPM	4
Rate	30
Tidal Volume	40 ml
W/Reservoir	100%
W/O Reservoir	85%
Measured Volume	
Bag	240 ml
Reservoir	600 ml

Latex-Free CPR Bags including Infant, Child, Small Adult and Adult

Mercury Medical's Latex-Free CPR Bags

6/Box

Adult CPR Bags



ORDERING INFORMATION

Product Number	Cush. Mask	Cont. Mask	O ₂ Resv.	Aero Resv.	Peep Valve	Pop-Off	Lock Clip	Omni-Link	Other CPR Parts	Manometer	Exhalation Filter
10-55000*		X	X								
10-55002*											
10-55003			X								
10-55004	X		X								
10-55007	X			X							
10-55008*				X							
10-55023*	X		X						1,2		
10-55024*	X			X	X						
10-55025*	X		X		X						
10-55027*	X		X					X			
10-55028*			X					X			
10-55045*	X		X			X	X		16		
10-55048*	X		X								X
10-55059*	X		X		X					X	
10-55061*	X			X				26			X
10-55068	X		X								
<i>ADULT SIZE WITH EXPANDABLE AEROSOL RESERVOIR</i>										<i>12/Box</i>	
10-55700									23,24		
10-55701	X								23,24		
10-55702					X				23,24		
10-55703*	X				X				23,24		
10-55704*	X					X	X		23,24		
10-55705*	X								23,24		X
10-55713*	X								23,24	X	

*Custom configurations require 10 business days to manufacture and a minimum order of 10 boxes per order.

Adult CPR Bags
with StatCO₂[®]
End Tidal CO₂ Detectors



Product Number	Cush. Mask	O ₂ Resv.	Aero Resv.	Expandable Aerosol Reservoir	StatCO ₂ End Tidal CO ₂ Detector	6/Box	12/Box
10-55064		X			X	X	
10-55065	X	X			X	X	
10-55066			X		X	X	
10-55067	X		X		X	X	
10-55711				X	X		X
10-55712	X			X	X		X

Mercury Medical's Latex-Free CPR Bags

6/Box
Small Adult CPR Bag



ORDERING INFORMATION

Product Number	Small Adult Cush. Mask	Cont. Mask	O ₂ Resv.	Aero. Resv.	Peep Valve	Pop-Off	Lock Clip	Omni-Link	Other CPR Parts	Manometer
10-55901	X								23	
10-55907	X			X						
10-55904	X		X							
10-55909						X	X		2, 23	
10-55910			X						2	X

6/Box
Child CPR Bag



ORDERING INFORMATION

Product Number	Cush. Mask	Cont. Mask	O ₂ Resv.	Aero. Resv.	Peep Valve	Pop-Off	Lock Clip	Omni-Link	Other CPR Parts	Manometer
10-55100*		X	X			X	X			
10-55102*						X	X			
10-55103			X			X	X			
10-55104	X		X			X	X			
10-55106				X		X	X			
10-55107	X			X		X	X			
10-55113*	X		X		X	X	X			
10-55115*				X		X	X		3	
10-55119*	X			X	X	X	X			
10-55122*	X		X			X	X		3	
10-55142*	X		X			X	X		3,5	
10-55150*	X		X		X	X	X		3	
10-55156*						X			3	
10-55163*	X		X			X			3	
10-55166*	X		X			X	X		1	
10-55167*	X		X			X				
10-55168*	X		X			X	X		3	X
10-55170*	X		X			X	X		3,5	X
10-55171*	X		X		X	X	X			X
10-55172*	X			X						
10-55173*	X		X			X	X		18	
10-55175*	X			X	X	X	X		3	X
10-55177*	X			X		X	X			X
10-55179*	X			X		X	X		3,5	X
10-55180*			X			X	X		3	X
10-55182*				X		X	X			
10-55183*	X		X			X	X		26	
10-55186*	X		X			X	X			X
10-55187*	X		X			X	X		26	
10-55188*				X		X	X		3	X
10-55189*	X			X		X	X		18	X
10-55190*			X			X	X			X

Mercury Medical's Latex-Free CPR Bags

*StatCO₂ End Tidal CO₂ Detector Foil Pack added to CPR drawstring bag.
The StatCO₂ has a 2-year expiration date
from date of manufacture.*

*Custom Configurations require 10 business days to manufacture and a minimum order of 10 boxes per order.

Child CPR Bag
with StatCO₂[®]
End Tidal CO₂ Detectors



ORDERING INFORMATION

Product Number	Cush. Mask	O ₂ Resv.	Aero Resv.	StatCO ₂ End Tidal CO ₂ Detector	6/Box
10-55185	X	X		X	X



6/Box
Infant CPR Bag

Product Number	Cush. Mask	Cont. Mask	O ₂ Resv.	Aero Resv.	Peep Valve	Pop-Off	Lock Clip	Omni-Link	Other CPR Parts	Manometer	Mini StatCO ₂
10-55200*		X	X		X	X					
10-55203			X			X	X				
10-55204	X		X			X	X				
10-55217*	X		X			X	X		2,5		
10-55223				X		X	X				
10-55224				X		X	X		5,20		
10-55229	X			X		X	X				
10-55234	X		X			X	X		5		
10-55241*	X		X		X	X	X				
10-55249*	X			X		X	X		5		
10-55250*	X			X		X	X		5,13		
10-55253	X			X		X			5,13		
10-55257*	X		X			X	X		5,8,21		
10-55264	X		X			X			5		
10-55266*						X	X		Size, 01 mask		
10-55271*			X			X					
10-55278*	X			X		X	X		5,21		
10-55280*	X			X		X	X		5	X	
10-55282*	X			X		X	X		11	X	
10-55283*	X			X		X	X			X	
10-55284*	X		X			X	X			X	
10-55287*			X			X	X			X	
10-55288*	X		X		X	X	X		5		
10-55289*			X			X			25		
10-55290*				X		X	X		Size 00/01 Masks		
10-55291*				X		X	X			X	
10-55292*	X		X		X	X	X			X	
10-55293*	X			X		X				X	
10-55294*			X			X	X		Size 00/01 Masks	X	
10-55295*	X		X			X			5,25		
10-55296*	X		X			X	X		18	X	
10-55297*	X		X		X	X	X		5	X	
10-55298*	X		X			X				X	
10-55299*	X		X			X	X		5	X	
10-552100*	X			X	X	X	X		5	X	
10-552101*	X			X	X	X	X			X	
10-552102*	X		X			X	X		26		
10-552107*				X		X	X			X	
10-552112*	X		X			X	X		5	X	X
10-552113*				X	X	X	X		27		

Mercury Medical's Latex-Free CPR Bags

OTHER CPR PARTS

- 1 Small Adult Mask
- 2 Child Mask
- 3 Infant Cushion Mask
- 5 Neonate Mask
- 6 Infant Patient Valve
- 7 Glued Oxygen Tubing
- 8 Glued Reservoir
- 9 Pediatric Mask
- 10 Standard Oxygen Tubing

- 11 Assembled except Aerosol Tubing
- 13 Removable Aerosol
- 15 Strapped O₂ Tube to Aerosol
- 16 Disconnect Wedge
- 18 Exhalation Filter
- 19 Bag Hanger
- 20 Hook Ring
- 21 Airway Pressure Line
- 23 Expandable Aerosol Reservoir
- 24 Packaged 12/case
- 25 Assembled except Oxygen Tubing
- 26 7' O₂ Tubing
- 27 Custom Zip Lock Bag with Handle

- All CPR Bag Patient Valves listed have built-in PEEP Adapters.
- All CPR Bag configurations come fully assembled.

CPR ACCESSORIES

- | | |
|----------|---------------------------------|
| 10-55310 | PEEP Adapter |
| 10-55330 | PEEP Valve |
| 10-55315 | PEEP Adapter & Valve Package |
| 10-55320 | Manometer Port |
| 10-55340 | Disconnect Wedge |
| 10-55360 | Omni-Link 12" |
| 10-55361 | Omni-Link II 6" |
| 10-55350 | Manometer (Single-Patient Use) |
| 10-55351 | Manometer with Tubing and Stand |
| 10-55375 | Exhalation Filter |

KEY TO TERMS USED

- | | |
|--------------------|---------------------------------|
| O2 Resv. | Oxygen Bag Reservoir |
| Cush. Mask | Cushion Mask |
| Cont. Mask | Contour Mask |
| Aero. Resv. | Aerosol Tubing Reservoir |